

# Whole Milk Powder

## DIGESTION OF WHOLE MILK POWDER USING COLDBLOCK™ DIGESTION TECHNOLOGY

*Author: Lorie-Anne Doig*

### Introduction

This application note will focus on the digestion of Whole milk powder using ColdBlock™ Digestion CB12L Technology.

### Method

A sample of Whole milk Powder was digested using the following method:

- 1.0g of sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 10mL H<sub>2</sub>SO<sub>4</sub> was added.
- Sample was digested at 100% power for 10 minutes.
- Chiller temperature was set to -5°C.
- After 10 minutes of digestion, 3-5 mL of H<sub>2</sub>O<sub>2</sub> (30%) was added until the solution was clear, and the reaction stopped.
- Sample was allowed to cool.
- Sample was normalized to 50mL using 2% HNO<sub>3</sub>.

### Instrument

ColdBlock™ Digestion CB12L Technology.

### General

This procedure is specific for the sample digested and may need modification for different samples to achieve the desired result.